# TITAN 9500 LED Headlight System

## Operation Manual SSL-9500

Sunoptic Technologies®  
6018 Bowdendale Avenue  
Jacksonville, FL 32216 USA

Customer Service: 904 737 7611  
Toll Free 877 677 2832

AJW Technology Consulting GmbH  
Breite Straße 3  
40213 Düsseldorf, Germany  
Telephone: [+49 211 54059 6030](tel:+49211540596030)

## TABLE OF CONTENTS

1. INTENDED USE
2. GENERAL WARNINGS
3. ASSEMBLY
4. MAINTENANCE
5. CLEANING
6. DESKTOP SMART CHARGER AND CALIBRATOR
7. REPLACEMENT PARTS
8. USER SERVICE
9. WARRANTY AND REPAIR
10. END OF PRODUCT LIFE
11. CHART OF MEDICAL DEVICE SYMBOLS USED

## 1. INTENDED USE

This LED TITAN 9500 Headlight system is designed to deliver illumination from a high intensity LED for surgical site illumination.

### **Congratulations on the purchase of your new LED Headlight system.**

The LED Headlight TITAN 9500 is a battery powered portable headlight with an adjustable spot and uses encrypted battery technology that requires the use of approved batteries obtained from the manufacturer or authorized distributors.

Additionally, a low battery alarm is offered in this system and will beep when power is low. Decreasing the output may extend battery life and therefore the amount of warning time before power down. This alarm is a convenience and the user should plan for backup lighting for any application using this system; see General Warnings.

## 2. GENERAL WARNINGS

**Note:** You must charge the batteries fully before using the headlight for the first time.

- The user of this product should be thoroughly familiar with the use and care of this product.
- The user should carefully study the manual before making any attempt to use this product clinically.
- Before each procedure, carefully check the battery state of charge to assure optimal duration for this procedure.
- LED high intensity headlights will produce some heat in addition to light. The high output of these headlights may become warm to the touch. This is normal operation.
- Light source can cause damage to the eye if viewed directly with unprotected eye.
- User is responsible for determining if interruption of light output will create an unacceptable risk. If this determination is made, arrangements (such as standby illuminator) should be made to reduce risk.
- Users are responsible for providing backup lighting systems for your application when using this device.
- The illuminator produces a highly concentrated (luminous energy per unit area) light beam. Qualified personnel must determine a safe working distance and intensity setting for each application.
- Not suitable for use in presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide.
- To prevent fire or electric shock, do not open or expose the illuminator unit to rain or moisture. Refer all servicing to qualified personnel only.
- The headlight should only be used with approved batteries supplied by Sunoptic Technologies. To order additional batteries, please call customer service 1-877-677-2832
- User should ensure that the air vents located on the headlight module are not obstructed or blocked to allow the headlight to receive the necessary cooling to prevent an over temperature shut off.
- This device meets CISPR 11 Class A limits and is suitable for use in a hospital environment. Performance of this device may be affected if in proximity of another device and/or equipment capable of producing high levels of RF transmissions. In the event performance of this device is affected due to high levels of RF transmissions, relocation of the suspected device and/or equipment producing high levels of RF transmissions

may eliminate the problem. If relocation of the suspected device and/or equipment producing the high level of RF transmission cannot be performed, relocation of the headlight holster position may reduce or eliminate the problem.

**NOTE:** This device may interfere with other electrical equipment if used outside a hospital setting.

- The headlight is permitted for use in ambient temperatures ranging from 0 to 30°C.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### 3. ASSEMBLING THE HEADLIGHT SYSTEM

*Before continuing remove the protective lens cap.*

Insert a fully charged encrypted battery into the Battery Holster and ensure it is fully seated down on the connector at the bottom of the battery cavity.

Should a non-encrypted battery be installed, two beeps will be heard upon insertion indicating that an incorrect battery is being used and the unit will not power up. Insert an encrypted battery and the unit will power up.

To attach the power cable from the LED Module to the Battery Holster, insert the push-pull quick connector of the LED Module into the receptacle on the lead from the Battery Holster until it clicks in place as shown here. This connection will not come off when positively engaged. You can rotate the connectors to find the keyway and then press together and not damage the pins.

The On/Off switch is also the LED intensity control. Simply rotate the knob clockwise or counter-clockwise to produce the desired light output.

The headlight has an iris to allow adjustment of the beam spot diameter. Rotate the fluted iris adjustment ring on the headlight clockwise or counterclockwise for the appropriate spot diameter.

The battery holster is worn on a belt; on the rear side of the holster is a metal belt clip that slips over and secures to a belt. The power cable secures into the cable clips on the side and rear of the headband.

To remove the LED Module from the Battery Holster, press the small raised button on the quick disconnect device located between the LED Module and the Battery Holster and pull straight out without twisting.


Next, place the headband on your head; you can make size adjustments by turning the knob on the rear of the headband until it feels comfortable. For more comfort you can adjust the cross band.

Once the headlight module and cable are connected, the unit is ready for use.

The system comes with gown clips to relieve the weight of the cable from the headband. One gown clip is attached below the shoulder.

Once the battery is depleted into the last shown “block” on the power gauge, the unit may beep once to indicate low battery levels. The unit may beep twice to indicate further reductions in power before shut-down.

 **CAUTION:** To avoid over-discharging of the battery pack, place the battery pack in the supplied charger once the holster begins to notify of low charge. **DO NOT** store batteries with low charge for extended periods of time, and **DO NOT** leave low charge batteries in the holster.

#### 4. MAINTENANCE

Storing the headlight system in a safe place will prolong the life of the unit. The following guidelines will help in sustaining the headlight system useful life:

- Keep the optical faces of the LED Headlight lens from touching hard surfaces, which may cause scratches on the surface. Scratches on the surface will diminish the light output. Using the lens cap during transport will help protect the lens.

#### 5. CLEANING

The headlight module and headband can be wiped down with commercially available cleansers commonly used for disinfection of electronic equipment in hospitals such as ethyl or isopropyl alcohols, disinfecting sprays containing quaternary ammonium compounds, or hydrogen peroxide. The lens should be cleaned only with lens tissue, available in any camera store; follow the directions on the package.

Do not use strongly caustic or acidic cleansers such as “Clorox” hypochlorite bleach, ammonia, muriatic acid or similar products. Do not use acetone, methyl ethyl ketone or halogenated / chlorinated hydrocarbon solvents or cleansers containing any of these restricted compounds.

Apply cleaning agents by light spray or dampened towels. Do not pour liquids onto the device. Do not allow liquids to enter the device seams or ventilation openings.

Follow all applicable bloodborne pathogen procedures as required by OSHA and/or your hospital when cleaning and disinfecting the product.

**Do not autoclave the LED Headlight module, headband, or battery.**

## 6. DESKTOP SMART CHARGER AND CALIBRATOR

Please refer to the separate battery and charger IFUs for detailed instructions on how to safely operate these devices.

## 7. REPLACEMENT PARTS

- For replacement parts, additional or optional equipment: Please call Customer Service @ 877 677-2832

## 8. USER SERVICE

The LED Headlight has no user or field serviceable components. It can only be serviced at the factory by factory trained technicians.

## 9. LIMITED WARRANTY

Your LED Headlight and battery holster carry a three year warranty from the date of shipment on workmanship and all defects of material.

Should your product prove to have such defects **Sunoptic Technologies®** will repair or replace the product or component part without charge. Please contact **Sunoptic Technologies®** for return authorization documentation.

You should carefully pack unit in a sturdy carton and ship it to the factory. Please include a note describing the defects, your name, telephone number and a return address. Warranty does not cover equipment subject to misuse, accidental damage, normal wear and tear or if transferred to a new owner without authorization from **Sunoptic Technologies®**. This warranty gives you specific legal rights and you may also have other rights that vary from state to state.

**POST WARRANTY REPAIRS:** You may return your product(s) for repair, shipping prepaid to the factory. Your product will be inspected and an estimate of repair charges will be submitted to you for approval. Payment must be received before repairs are completed.

- In the U.S call: 877 677 2832 (toll free)
- Customer Service 904 737-7611

## 10. END OF PRODUCT LIFE

In accordance with the European Waste from Electrical and Electronic Equipment (WEEE) directive, we encourage our customers to recycle this product whenever possible. Disposal of this unit must be performed in accordance with the applicable local environmental regulations.

In the US a list of recyclers in your area can be found at: <http://www.eiae.org/>.

Please contact customer service to issue a return authorization to return product to manufacturer at the end of product life.

## 11. CHART OF MEDICAL DEVICE SYMBOLS USED

|   |   |
|---|---|
|    | CE mark   |
|    | Manufacturer  |
|    | "Authorized Representative in the European Community."  |
|    | Caution, consult accompanying documents                 |
|    | Do not dispose of in ordinary municipal waste, Recycle. |
|    | Medical device  |
|    | Non-Sterile   |
|    | Storage / Shipping Temperature                          |
|    | Storage / Shipping Humidity                             |
|   | Barometric Pressure                                     |
|  | Unique Device Identifier                                |
|  | Keep Dry  |

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