[illegible]

CUTANPLAST[®]

Product description

Cutanplast[®] is a sterile gelatine sponge from animal origin, insoluble in water, easy to handle and re-absorbable. Its purpose is to facilitate the stopping of bleeding. The sponge is whitish with uniform porosity.

Action mechanism

Once saturated with blood, the product promotes coagulation by serving as a coadjutant for controlling local haemorrhage. If utilised correctly and in suitable amounts, Cutanplast[®] is reabsorbed within 3-4 weeks. The amount of product to be used and its re-absorption depend on the extension of the bleeding surface and the type of tissue on which it is applied.

Indications

Cutanplast[®] can be used dry or soaked with sterile saline solution, for surgical interventions, except for ophthalmic ones, through simple application or by exerting slight pressure like a haemostatic device for controlling the bleeding of venous or arteriolar capillaries, and in cases where other conventional haemostatic procedures (stitching, compression, binding etc.) cannot be carried out.

Cutanplast[®] can be cut to the desired size without crumbling, using sterile scissors, and it has low adhesiveness with respect to forces and surgical instruments.

Using Cutanplast[®] Powder requires mixing with sterile saline solution.

Contraindications

Do not use Cutanplast[®] for closing skin cuts as it could interfere with healing and cicatrization.

Do not use Cutanplast[®] in intravascular cavities due to the risk of embolization.

Do not use Cutanplast[®] for embolization procedures.

Do not use Cutanplast[®] in patients with ascertained allergy towards the product.

Do not use Cutanplast[®] on infected wounds.

Beware

- The use of Cutanplast[®] does not authorise abandonment of normal haemostatic intervention.
- After attaining haemostasis, if possible remove the product from the site of application to avoid the risk of displacing the product and possible compression of the adjacent structures
- Remove Cutanplast[®] from the site of application when used on bone holes, in the spinal cord and/or in the nerve or optical chiasm or in the areas surrounding or adjacent thereto.
- Remove Cutanplast[®] from the site of application if used in neurosurgery or in the areas surrounding or adjacent to the nerves given that by increasing the volume the product may cause compression leading to damage of nervous tissue.
- Remove Cutanplast[®] in cases of strong arterial haemorrhage, in cases of blood stagnation or stagnation of other organic liquids where the source of haemorrhage is entirely occluded.
- Do not use Cutanplast[®] in post-partum interventions and for controlling haemorrhage.
- The safety and efficiency of Cutanplast[®] in ophthalmic interventions is not verified
- The safety and efficiency of Cutanplast[®] on children and expectant women is not verified.

Precaution measures

Beware: Cutanplast[®] is supplied sterile and disposable. Unutilised fractions of the product should be discarded. Any opened and unutilised packaging should be discarded.

Warning: When in use on closed or narrow tissue spaces preliminary compression of the product is advisable to prevent its increase volume-wise, due to contact with organic fluids, causing damage.

Warning: When indicated to fill a cavity to obtain haemostasis, Cutanplast[®] shall not be used this way unless all surplus product is removed after achieving haemostasis.

Warning: Use Cutanplast[®] in minimal amounts to reach haemostasis and always remove any surplus product from the site of application.

Warning: Do not use Cutanplast[®] together with autologous blood recovery systems, given that it has been proven that particles of haemostatic agents could pass through the transfusion filters of the purification systems.

Warning: Do not use Cutanplast[®] together with methyl methacrylate cement. It has been proven that the porous structure weakens the cement used for fixing the prosthesis to the bone surface.

Warning: There have been cases of incomplete re-absorption and loss of heating due to the use of haemostatic sponges in eardrum-plastic operations.

Warning: Do not use Cutanplast[®] in cases of coagulation disorders.

Warning: The association of Cutanplast[®] with pharmacologically active substances or other substances shall be evaluated by the surgeon and the responsibility thereof.

Warning: The safety and efficiency of Cutanplast[®] in urology surgical procedures is not verified. Cutanplast[®] should not be left in renal pelvis, in renal calyces, in the urethra or in the ureter so as to avoid the possible formation of calcifi.

Warning: CUTANPLAST[®] is not recommended for treating cases of epistaxis given that the product could be inhaled thus causing respiratory problems.

How to use

Carefully read the present use instructions for use before using the product

Ensure that the packaging is whole before using. Sterility is guaranteed only if the packaging is not damaged or wet.

Be very cautious when it comes to aspects where handling is particularly important.

Always use the minimal amount of product to attain haemostasis.

Cutanplast[®] can be used dry as it is soaked with saline solution.

Apply the product on the area subject to haemorrhage by applying slight pressure until haemostasis is attained.

If the product cannot be removed upon attaining haemostasis, carefully eliminate the surplus used product.

Use the product by the expiry date indicated on the primary packaging, secondary packaging and on the box.

The product is supplied as per the formats indicated in the present instructions leaflet.

Cutanplast[®] Powder:

Pour the powder into a suitable sterile container and add 3-6 ml of sterile saline solution.

Dispensing of the powder is attained by compression using fingers, suitably covered with sterile gloves, on the bottom of the container up to obtain a malleable paste.

The paste thus obtained may be spread, rolled or applied according to preference.

Preservation

Cutanplast[®] shall be stored at room temperature and in a dry place

Cutanplast[®] shall be stored away from direct sunlight

Cutanplast [®] Standard	70 x 50 x 10 mm	Ref. 05170502
Cutanplast [®] Special	70 x 50 x 1 mm	Ref. 05380502
Cutanplast [®] Anal	80 x 30 mm	Ref. 05480502
Cutanplast [®] Dental	10 x 10 x 10 mm	Ref. 05510101
Cutanplast [®] Film	200 x 70 x 0,5 mm	Ref. 05770202
Cutanplast [®] Film x 50	200 x 70 x 0,5 mm	Ref. 05770203
Cutanplast [®] Dial	30 x 30 x 10 mm	Ref. 05230302
Cutanplast [®] Small	50 x 30 x 10 mm	Ref. 05680302
Cutanplast [®] Large	80 x 125 x 10 mm	Ref. 058801251
Cutanplast [®] 20 x 60	20 x 60 x 7 mm	Ref. 05520602
Cutanplast [®] 80 x 20	20 x 80 x 10 mm	Ref. 05580202
Cutanplast [®] 80 x 50	80 x 50 x 10 mm	Ref. 05190502
Cutanplast [®] Powder		

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