OKCEL®

Oxidized Resorbable Cellulose

HAEMOSTAT

Synthesia, a.s.

Semtín 103
530 02 Pardubice
Czech Republic

www.synthesia.eu
www.okcel.eu

Návod k použití
Gebrauchsanweisung
Instructions for use
Indicaciones de uso
Notice d'instruction
Upute za upotrebu
Istruzioni per l'uso
Használati Utasítás
Instrukcja używania
Instrucțiuni de utilizare
Οδηγίες Χρήσης
Návod na použitie
Uputstvo za upotrebu
Bruksanvisning
Kullanma talimatları
Правила применения
Інструкції з використання
تعليمات الاستخدام

OXY
cellulose

18640_I_O_G_2

CE 1023

Čeština	1
Deutsch	5
English	9
Español	13
Français	17
Hrvatski	21
Italiano	25
Magyar	29
Polski	33
Română	37
Ελληνική	41
Slovenčina	45
Srbski	49
Svenska	53
Türkçe	57
Русский	61
Українски	65
69	عربي

A	B	C	D
OKCEL® H-T	H-T 101	1.5 x 1.5 cm	10 units of 10 pcs
	H-T 151	1.5 x 1.5 cm	20 units of 2 pcs
	H-T 501	5 x 1.25 cm	15 pcs
	H-T 507	5 x 7 cm	15 pcs
	H-T 510	7 x 10 cm	15 pcs
	H-T 535	5 x 35 cm	10 pcs
	H-T 540	10 x 20 cm	10 pcs
OKCEL® H-D	H-D 202	2.5 x 2.5 cm	15 pcs
	H-D 209	2.5 x 9 cm	15 pcs
	H-D 575	5 x 7.5 cm	10 pcs
	H-D 710	7 x 10 cm	10 pcs
	H-D 1420	14 x 20 cm	10 pcs
	H-D 1523	15 x 23 cm	10 pcs
OKCEL® F	F 205	2.5 x 5 cm	10 pcs
	F 575	5 x 7.5 cm	10 pcs
	F 510	5 x 10 cm	10 pcs
	F 1010	10 x 10 cm	10 pcs
	F 1020	10 x 20 cm	10 pcs
OKCEL® S	S 205	2.5 x 5 cm	10 pcs
	S 505	5 x 5 cm	10 pcs
	S 510	5 x 10 cm	10 pcs
	S 1010	10 x 10 cm	10 pcs

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OKCEL® Oxidized Resorbable Cellulose

Description:

OKCEL® is a sterile absorbable knitted fabric (OKCEL® H-T and OKCEL® H-D) or a sterile absorbable soft layered material (OKCEL® F and OKCEL® S), prepared through the controlled oxidation of cellulose. The product's colour is white to slightly yellowish and it can be cut without fraying out.

The OKCEL® H-T product line with regular knitted fabric density is used to control capillary, minor venous and minor arteriolar diffuse bleeding.

The higher density knitted fabric OKCEL® H-D providing greater thickness and improved efficiency is used to control haemostasis in higher volume capillary and venous or arteriolar bleeding.

The non-woven fibrous version OKCEL® F with extremely high flexibility is used to control haemostasis over a large area, for surface applications on irregularly shaped bleeding sites or in areas that are difficult to access. It can be easily placed on the bleeding site using forceps. Any amount of OKCEL® F can be easily separated in such a way that the fibres continue to hold together without uncontrolled release into the operative site.

The non-woven reinforced fibrous version OKCEL® S with reduced weight and increased strength is used to control capillary, venous and minor arteriolar diffuse bleeding and it may be more suitable for endoscopic use. The material maintains its original structure even when saturated with blood and can be easily manipulated in the bleeding site with no signs of disintegration.

Mechanism of Action:

When saturated with blood, OKCEL® swells into a gelatinous mass, still maintaining its original structure. The product then helps in the formation of a clot by initial denaturation of blood proteins, thus supporting local haemostasis and haemorrhage control. When used properly in minimal amounts, OKCEL® is absorbed with practically no tissue reaction, depending on the degree of saturation with blood, and the tissue bed.

In addition to topical haemostasis, the bacteriostatic and bactericidal properties of OKCEL® are another positive effect, which have been proved in vitro. Low pH inhibits the growth and multiplication of both gram positive and gram negative microorganisms including aerobic and anaerobic bacteria. Nevertheless, OKCEL® is not intended to be a substitute for systemically administered therapeutic or prophylactic anti-microbial agents to control or prevent post-operative infections.

Indications:

OKCEL® is used adjunctively during surgical and minimally invasive procedures to control of haemostasis of capillary, venous and minor arteriolar bleeding, when conventional haemostatic methods such as ligation are not effective or practical. OKCEL® can be cut to size for use (not only) in endoscopic procedures.

Available Forms:

OKCEL® H-T (regular density knitted fabric) is supplied sterile in a double pack.

OKCEL® H-D (high density knitted fabric) is supplied sterile in a double pack.

OKCEL® F (regular density fibrous version) is supplied sterile in a double pack.

OKCEL® S (reduced weight reinforced fibrous version) is supplied sterile in a double pack.

The imprint of the Instructions for Use contains a table with manufactured dimensions and configurations of the sales packaging of each of OKCEL® forms as follows: column A – medical device form, column B – catalogue number, column C – dimensions, column D – number of units a number of pieces in the unit, event. number of pieces (pcs).

Instructions for Use:

Check the integrity of the packaging before using OKCEL®. Maintain product sterility when removing OKCEL® from the packaging and handling it.

The outer pouch should be peeled open to allow the inner pouch to drop on a sterile tray or towel. The inner pouch should be peeled open by a gloved member of the operating team under sterile conditions. The product can then be removed using forceps. The required amount of OKCEL[®], adjusted as necessary with scissors, is laid on the bleeding site or held firmly against the tissue until bleeding stops. This haemostatic agent is easy to handle, therefore the position of OKCEL[®] can be altered if necessary. Any open unused product should be discarded. Before surgical closure remove any excess amount to facilitate biological dissolution and minimise a possible reaction to the foreign body.

Contraindications:

OKCEL[®] must not be used for implantation into bone defects, such as fractures, because of the potential of interfering with callus formation and a theoretical chance of cyst development.

OKCEL[®] must not be used in conjunction with methyl methacrylate-based tissue adhesives, for example in orthopaedic surgery, because their presence may reduce the bonding agent's adhesion to the bone.

OKCEL[®] must not be used for implantation around the optic nerve and the chiasma opticum and must not remain in contact with other neural pathways, where increased pressure induced by the swelling of the absorbable haemostatic material may affect the neural pathway function.

OKCEL[®] must not be used to control bleeding from large blood vessels.

OKCEL[®] must not be used as prevention to avoid adhesions.

OKCEL[®] must not be used on non-haemorrhagic serous oozing surfaces, since body fluids other than whole blood, such as serum, do not react with the product to produce satisfactory haemostatic effect.

OKCEL[®] do not substitute careful surgical procedures, such as sutures and ligatures.

Enclosing the haemostat in contaminated undrained wound may lead to complications and such situations must be avoided.

OKCEL[®] must not be used inside the blood circulatory system. Circulating fragments of the haemostat may lead to partial or complete obstruction of the vascular system.

When used in vascular surgery, extra care must be taken to avoid penetration of the haemostat into the

vascular system. Otherwise, loose fragments of the product may lead to partial or complete obstruction of the circulatory system.

OKCEL[®] is not designed, supplied or intended for any use other than as indicated. It is the physician's responsibility to determine whether any physical impairment of the patient could contraindicate the use of this product.

Warnings:

OKCEL[®] may be used by trained professionals only. OKCEL[®] is intended for single use only. Do not reuse or resterilize! Discard any open unused material. Do not use if the product is yellow/brown. Brownish colour of the haemostat and a tendency to crumble may indicate damage to the product caused by incorrect storage conditions.

Do not use if any of the double pack layers or their content are damaged or appear to be damaged.

Wadding or tight packing with topical absorbable haemostats must be avoided, especially within the bone compartment of the central nervous system or within other rigid cavities, where increased pressure induced by the enlarged volume of the absorbable haemostatic material may affect function of nerve tissue and may even cause its necrosis.

If used as a wrap in vascular surgery, do not apply too tight.

When OKCEL[®] is used to achieve haemostasis within, around, or in the proximity of foramina in bones, in areas confined by bones, or in the proximity of the optic nerve and the chiasma opticum, the haemostat must always be removed after haemostasis has been achieved; otherwise the material will swell, exert unwanted pressure and compress the tissue.

Regardless of the type of surgical procedure, doctors must consider carefully whether it is appropriate to remove the haemostatic agent once haemostasis is achieved.

OKCEL[®] is to be applied dry and must not be moistened prior to application! Make sure to remove as much blood as possible from the bleeding site before applying a local absorbable haemostat.

OKCEL[®] must not be impregnated with anti-infective agents or with other materials such as buffer solutions or other haemostatic substances. Haemostatic effect is not enhanced by the addition of thrombin. Thrombin activity is suppressed by the low pH of oxycellulose. In general, when Redon drainage is used,

it must not be placed close to the applied haemostat in order to ensure safe drainage and prevent plugging. Doctors should be aware that cases have been reported of the formation of encapsulated fluid collections and reactions to a foreign body in this type of product. This should be taken into consideration during further monitoring or diagnostic imaging as a possible cause of certain clinical findings which could otherwise indicate a pathological condition.

Precautions:

Use only as much of the product as necessary for achieving haemostasis and for holding it firmly in place until permanent haemostasis is achieved.

When used in urological cases, extra care must be exercised to prevent the device from plugging the urethra, ureters or a catheter.

OKCEL® must not be used as a surface dressing as it inhibits the growth of new skin.

If OKCEL® is used temporarily to treat bleeding in large open wounds, it should be placed so as not to overlap skin edges.

OKCEL® must not be used in infected areas.

The use of OKCEL® should not be preceded by the application of caustic chemicals, such as silver nitrate or other similar products, since the absorption of the haemostat may be reduced in chemically cauterized areas.

As a passive haemostat OKCEL® is suitable for use in patients with intact coagulation system. There is no clinical data on the use of the device in patients under 12 months of age

Potential Complications:

Potential complications include, but are not limited to:

- Stenosis
- Paralysis, nerve damage and/or blindness
- Drainage complications after cholecystectomy and urethral procedures
- "Burning" and "stinging" sensations, sneezing
- Risks normally associated with interventional procedures

There have been reports of a stenotic effect when oxidized cellulose has been applied as a wrap during vascular surgery. Although it has not been established that the stenosis was directly related to the use of oxidized cellulose, it is important to be cau-

tious and to avoid applying the material tightly as a wrapping.

Paralysis and nerve damage have been reported when oxidized cellulose was used around, within, or in the proximity of foramina in bone, areas of bony confine, the spinal cord, and/or the optic nerve and chiasma opticum. While most of these reports have been in connection with laminectomy, reports of paralysis have also been received in connection with other procedures. Blindness has been reported in connection with surgical repair of a lacerated left frontal lobe when oxidized cellulose was placed in the anterior cranial fossa. Possible prolongation of drainage in cholecystectomy and difficulty in passing urine via the urethra after prostatectomy have been reported. A blocked ureter after kidney resection has been reported, which required postoperative catheterization. Occasional reports of "burning" and "stinging" sensations and sneezing when oxidized cellulose has been used as packing in epistaxis are believed to be due to the low pH of the product.

A burning sensation has been reported when oxidized cellulose was applied after nasal polyp removal and after haemorrhoidectomy. Headache, burning, stinging and sneezing in epistaxis and other rhinology procedures, as well as stinging where oxidized cellulose was applied on surface wounds (varicose ulcerations, dermal abrasions and donor sites) have also been reported.

Interactions:

Not known.

Overdose/intoxication:

Not known.

Effects on the ability to drive and operate machinery:

Not known.

Use in pregnancy and lactation:

Not known.

Storage and shelf life:

Keep dry. For storage temperature conditions please refer to the label on the box and pouches. Protect from direct sunlight. Keep in the original packaging. Do not use beyond the use by date stated on the label. Browning of the product and a tendency to

crumble indicate incorrect storage. Do not use the product if damaged! Sterilized by gamma irradiation. Do not resterilize.

Recommendation for disposal:

Dispose of medical device by authorized person.

1. Transport packaging and sale packaging, including accompanying documents (if not required) – sorted waste (15 01 01¹ paper and cardboard packaging).
2. Medical device including primary and secondary packaging, residues of medical devices used and unused, which were handled in the operating room – infectious waste (18 01 03 wastes whose collection and disposal is subject to special requirements in order to prevent infection).

3. Unused medical device, not handled in operating room (e.g. after use-by date)

- a) Transport packaging and sale packaging including accompanying documents – sorted waste (15 01 01 paper and cardboard packaging),
- b) primary and secondary packaging including medical device – 18 01 04 wastes whose collection and disposal is not subject to special requirements in order to prevent infection.

Any serious incident related to this medical device should be reported to manufacturer and the appropriate competent authority.

Synthesia, a.s.
Semtín 103
530 02 Pardubice, Czech Republic

1023

REF

CATALOGUE NUMBER

LOT

BATCH CODE

SIZE

PRODUCT SIZE

STERILE R

STERILIZED USING IRRADIATION

USE BY DATE

CONSULT INSTRUCTIONS FOR USE

DO NOT RESTERILIZE

DO NOT REUSE

TEMPERATURE LIMITATION (STORAGE REQUIREMENTS)

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