



Xe

Replacement Lamps and Modules

CROSS REFERENCE LIST



MEDICAL OEM EQUIPMENT	MODEL NUMBER / DESCRIPTION	LAMP PART NUMBER	LAMP MODULE PART NUMBER
Applied Fiber Optics	Sunbeam Xenon Light Source	SSX1579	
Excelitas / Perkin Elmer / LLC	PE175BF	SSX175BF	SSX1339
Excelitas / Perkin Elmer / LLC	Cermax 300W, PE300BF	SSX300BF	SSX1338
BFW	Maxenon Xi300	SSX300BF	
Cooper Surgical	Endomed LSS-700	SSX300BF	
Dyonics (Smith & Nephew)	Dyo-Bright 3000XE, 300XL	SSX300BF	SSX1407
Gulf Fiberoptics	Visualux LX 300	SSX300BF	SSX1338
Isolux	XSB1300 300W Illuminator, Millenium	SSX300BF	SSX1338
Leica	M720F60		SSX1585
Leica	OH5, OH4, OH3, OHS1, MS3, MS2, MS1, F40	SSX300BF	SSX1426
Linvatec / Conmed	8430, 8180, 8170	SSX300BF	SSX1426
Linvatec / Conmed	LS7700	SSX1579	SSX1752
Luxtec / Integra	LX300	SSX300BF	

MEDICAL OEM EQUIPMENT	MODEL NUMBER / DESCRIPTION	LAMP PART NUMBER	LAMP MODULE PART NUMBER
Luxtec / Integra	9300XDP, 9300XSP	SSX300BF	SSX1739
Luxtec / Integra	9175	SSX175BF	
Nextec	300W Xenon, 1300	SSX300BF	
Olympus	XLS (175W) Y-1176 or CLV-A (175W)	SSX175BF	
Olympus	CLVS, CLV-10, CLV-F10, CLV-S20, CLV-S40, CLV-160, ILX6300	SSX300BF	
Olympus	CLV-190/CLV-290 (MAJ-1817)	SSX1854	
Olympus	CLV-180 (MD-631)	SSX1846	
Pentax	EPK-i	SSX1576	
Pentax	EPM3500, EMP3000, 3300, LX3000	SSX300BF	
Pilling Weck	300W Xenon, 521393	SSX300BF	SSX1338
Products For Medicine	White Sun 300W	SSX300BF	SSX1338
QED	3010 300W Light Source	SSX300BF	SSX1338
Smith & Nephew	300 XL (Also listed under Dyonics)	SSX300BF	

MEDICAL OEM EQUIPMENT	MODEL NUMBER / DESCRIPTION	LAMP PART NUMBER	LAMP MODULE PART NUMBER
Snowden Pencer	Xenon L/S 300W	SSX300BF	
Storz (175 Watt Systems)	615B, 615C, 20131520, 20132120	SSX175BF	
Storz	300(20133020) 20133206 Lamp Only)	SSX300BF	
Storz (175 Watt Module)	20132025	SSX175BF	
Storz (175 Watt Lamp)	20132026	SSX175BF	
Storz (300 Watt Systems)	487, 630, 201330-20	SSX300BF	
Storz (300 Watt Modules)	20133625, 20133027	SSX300BF	
Storz (300 Watt Lamp)	20133028	SSX300BF	
Stryker X8000	Module Swap-send in old module with IButton)	SSX1950	
Stryker	X6000	SSX300BF	
Stryker	X7000	SSX1492 Use only with SSX1571 Assembly	SSX1571 (use original housing)
Stryker	X8000	SSX1495 Use only with SSX1567 Assembly	SSX1567 (use original housing)
Sunoptic Surgical	SolarMaxx 180		SSX0048

MEDICAL OEM EQUIPMENT	MODEL NUMBER / DESCRIPTION	LAMP PART NUMBER	LAMP MODULE PART NUMBER
Sunoptic Surgical	SolarMaxx 300		SSX0030
Sunoptic Surgical	TITAN 180 Light Source		SSX0049
Sunoptic Surgical	TITAN 300 Light Source		SSX0035
Sunoptic Surgical	TITAN 400 Light Source		SSX0040
Sunoptic Surgical	TITAN X300 Light Source		SSX0300
Sunoptic Surgical	TITAN X350 Light Source		SSX0350
Sunoptic Surgical	TITAN X400 Light Source		SSX0400
Sunoptic Surgical	TITAN X450 Light Source		SSX0450
Sunoptic Surgical	Pro300 Light Source		SSXP300
Weck (Pilling)	FRX 2000, WECKLITE 300, 521393	SSX300BF	
Welch-Allyn	49820	SSX300BF	SSX1338
Welch-Allyn	Pro Xenon		SSX1681



Illuminate Excellence: Xenon Lamps for Precision Surgery

MEDICAL OEM EQUIPMENT	MODEL NUMBER / DESCRIPTION	LAMP PART NUMBER	LAMP MODULE PART NUMBER
Wolf	300W Xenon, 5131, 5151XE	SSX300BF	SSX1338
Zeiss	Superlux 300, Superlux 301 (with timer)	SSX300BF	SSX1426
Zeiss	Superlux 175 (with timer)	SSX175BF	

When it comes to sourcing replacement lamps and modules for surgical procedures, Sunoptic Surgical stands out as the premier choice. Proudly made in the USA, our products are crafted with the highest quality materials to ensure exceptional performance and reliability in the operating room. Our commitment to excellence is underscored by our ISO 13485 certification, guaranteeing adherence to the most stringent standards in medical device manufacturing. Depend on Sunoptic Surgical for unparalleled quality and service, enhancing the success and safety of your surgical procedures.

Designed, Developed, and Manufactured in the United States



www.SunopticSurgical.com

 904-737-7611



Customer.Service@SunopticTech.com

