

Reusable Surgical Instrument Cleaning Brushes



Exterior Brush
(Soft)

Order Code:
LT530



Exterior Brush
(Medium)

Order Code:
LT520



Exterior Brush
(Rigid)

Order Code:
LT510



Exterior Brush
(Aggressive)

Order Code:
LT510-C



Orthopaedic
Rasp Brush
(Standard)

Order Code:
LT550



Orthopaedic
Rasp Brush
(Aggressive)

Order Code:
LT550-C



Acetabular
Reamer
Cleaning Brush

Order Code:
LT540



Diathermy
Cleaning Brush
(Standard)

Order Code:
LT560



Diathermy
Cleaning Brush
(Aggressive)

Order Code:
LT560-C



Precision
Cleaning Brush
(Soft)

Order Code:
LT430



Precision
Cleaning Brush
(Rigid)

Order Code:
LT410



Precision
External
Cleaning Brush

Order Code:
LT570



Orthopaedic
Shaver
Cleaning Brush

Order Code:
LT565



Microflex
Channel
Cleaning Brush

Available in
Various Sizes



G-Flex
Refillable
Cleaning Brush


Available in
Various Sizes



Uniplex (UK) Ltd,
Sheffield, S3 7AF
Web: www.uniplexuk.com
Tel: +44 (0)114 241 3410
Email: info@uniplexuk.com



Single-Use Surgical Instrument Cleaning Brushes



Endoscopy Channel Cleaning Brushes

Device	Code:55cm	Code:75cm	Code:100cm	Code:150cm	Code:200cm	Code:250cm	Code:300cm	Colour
Channel Diameter: 1.2mm - 1.5mm	SU055-1215	SU075-1215	SU100-1215	SU150-1215	SU200-1215	SU250-1215	SU300-1215	Blue
Channel Diameter: 1.8mm - 2.4mm	SU055-1824	SU075-1824	SU100-1824	SU150-1824	SU200-1824	SU250-1824	-	White
Channel Diameter: 2.6mm - 3.2mm	SU055-2632	SU075-2632	SU100-2632	SU150-2632	SU200-2632	SU250-2632	-	Orange
Channel Diameter: 3.7mm - 4.8mm	SU055-3748	SU075-3748	SU100-3748	SU150-3748	SU200-3748	SU250-3748	SU300-3748	Red



Endoscope Port Cleaning Brush

Device	Code	Colour
Brush Diameter: 4mm & 8mm	SULV-08	Silver
Brush Diameter: 6mm & 11mm	SULV-11	Silver



Laparoscopic Channel Cleaning Brushes

Device	Code:55cm	Colour
Channel Diameter: 4mm - 5mm	SU-L05	Blue
Channel Diameter: 4.5mm - 5.5mm	SU-L55	White
Channel Diameter: 5.5mm - 7mm	SU-L07	Orange
Channel Diameter: 7mm - 9mm	SU-L09	Red
Channel Diameter: 10mm - 12mm	SU-L12	Blue
Channel Diameter: 12mm - 16mm	SU-L16	White



Channel Cleaning Brushes

Device	Code:30cm	Colour
Channel Diameter: 7.0mm	SU-M25	Blue
Channel Diameter: 9.0mm	SU-M30	White
Channel Diameter: 12.0mm	SU-M40	Orange
Channel Diameter: 15.0mm	SU-M50	Red



Flexible Endoscope Exterior Brush

Device	Code	Colour
Exterior Use	SU-BS	White



Tracheostomy Tube Cleaning Brushes

Device	Code:15cm	Colour
Channel Diameter: 6mm - 7mm	SU-CT01	Blue
Channel Diameter: 8mm - 9mm	SU-CT02	White
Channel Diameter: 10mm	SU-CT03	Orange
Channel Diameter: 11mm - 12mm	SU-CT04	Red



Uniplex (UK) Ltd,
Sheffield, S3 7AF
Web: www.uniplexuk.com
Tel: +44 (0)114 241 3410
Email: info@uniplexuk.com